



# Long-term Outcomes with Biodegradable Polymer Sirolimus-eluting Stents versus Durable Polymer Everolimus-eluting Stents in ST-segment Elevation Myocardial Infarction (STEMI): 5-year follow-up of the BIOSTEMI randomized trial

## Conclusions

- In patients with STEMI undergoing primary PCI, Orsiro is superior to Xience with respect to the rates of TLF at 5 years of follow-up, a difference driven by a numerically lower risk of clinically-indicated TLR.<sup>1</sup>
- Orsiro shows 31% significantly less target lesion failure at 5-year in STEMI patients: Orsiro DES: 7.7% vs. Xience DES: 11.1% [BIOSTEMI with historical information RR, 0,70; 95% BCI, 0.51-0.95, Bayesian posterior probability, 0.988]<sup>1</sup>

## Study design

Investigator-initiated, prospective, multicentre, assessor-blinded, randomized (1:1), controlled, superiority trial comparing Orsiro and Xience in STEMI patients undergoing primary PCI.

## Endpoints

### Primary Endpoint for BIOSTEMI

Target Lesion Failure (TLF) at 12-month follow-up defined as the composite of:

- Cardiac Death
- Target Vessel-Myocardial Re-Infarction (TV-reMI)
- Clinically Indicated-Target Lesion Revascularization (CI-TLR)

### Primary Endpoint for BIOSTEMI Extended Survival


TLF at 60-month follow-up defined as the composite of:

- Cardiac Death
- TV-reMI
- CI-TLR

### Selected Secondary Endpoints for both BIOSTEMI and BIOSTEMI ES

Individual components of the primary endpoint, All Cause Death, Target Vessel Revascularization (TVR), Target Vessel Failure (TVF), Definite Stent Thrombosis, Definite or Probable Stent Thrombosis (ST)

Patient characteristics <sup>1</sup>	Orsiro n = 649	Xience n = 651
Age, years*	62.2 ± 11.8	63.2 ± 11.8
Male	79%	73%
Active Smoker	46%	39%
Diabetes Mellitus	11%	13%
BMI [kg/m <sup>2</sup> ]*	26.9 ± 4.3	26.8 ± 4.3
Previous MI	4%	4%
Previous PCI	5%	5%
Previous CABG	0.3%	1%



## Angiographic and procedural characteristics<sup>1</sup>


	Orsiro n = 816**	Xience n = 806**
Number of lesions/patient*	1.26 ± 0.57	1.24 ± 0.52
Total Occlusion	49%	55%
Thrombus Aspiration	30%	31%
Baseline TIMI flow		
0 or 1	55%	59%
2	13%	14%
3	31%	27%
Cardiogenic shock	3%	3%
Small vessel (minimum stent diameter ≤3.0 mm)	36%	40%
Bifurcation treatment (including left main coronary artery)	12%	14%
Long Lesions (total stent length ≥20 mm)	71%	71%

\* Data shown as mean ± SD; \*\* Number of lesions


## Primary Endpoint – TLF at 5 years<sup>1</sup>

Orsiro – Proven safety and efficacy at short- and long-term follow-up with 31% significantly less TLF at 5-year in STEMI patients.<sup>3</sup>

TLF at 5 years<sup>1</sup>




TLF Landmark Analysis at 2 years<sup>2</sup>




## Selected Secondary Endpoints at 5 years

Orsiro is superior to Xience with respect to the rates of TLF at 5 years of follow-up, a difference driven by a numerically lower risk of clinically-indicated TLR.<sup>3</sup>

CI-TLR at 5 years<sup>1</sup>



CI-TLR Landmark Analysis at 2 years<sup>2</sup>



## Other secondary outcomes<sup>2</sup>

			BIOSTEMI with historical information from BIOSCIENCE		BIOSTEMI only Without historical information from BIOSCIENCE	
	Orsiro n = 649	Xience n = 651	Rate Ratio (95% BCI)	Bayesian Posterior Probability	Ratio Ratio (95% BCI)	Bayesian Posterior Probability
Target Lesion Failure	8%	11%	0.70 (0.51-0.95)	0.988	0.68 (0.47-0.98)	0.981
Cardiac Death	5%	6%	0.81 (0.54-1.23)	0.839	0.89 (0.55-1.43)	0.677
TV-reMI	2%	3%	0.76 (0.41-1.34)	0.833	0.67 (0.32-1.35)	0.868
CI-TLR	3%	5%	0.68 (0.40-1.06)	0.956	0.56 (0.32-0.96)	0.982
Target Vessel Failure	10%	13%	0.74 (0.55-0.97)	0.984	0.71 (0.51-0.98)	0.982
CI-TVR	4%	6%	0.59 (0.34-0.98)	0.979	0.56 (0.34-0.92)	0.990
POCE	16%	18%	0.88 (0.66-1.14)	0.836	0.87 (0.67-1.13)	0.847
Definite Stent Thrombosis	2%	3%	0.58 (0.28-1.18)	0.933	0.59 (0.28-1.20)	0.927

## Principal investigators

Dr. J.F. Iglesias, Geneva University Hospital, Geneva, Switzerland  
 Dr. O. Muller, Lausanne University Hospital, Lausanne, Switzerland  
 Dr. T. Pilgrim, Bern University Hospital, Bern, Switzerland

BCI: Bayesian credible interval, BPP: Bayesian Posterior Probability, CABG: Coronary Artery By-Pass Graft, CI: Confidence Interval, PCI: Percutaneous Coronary Intervention, RR: Risk Ratio.

1. Iglesias, JF. Long-Term Outcomes with Biodegradable Polymer Sirolimus-Eluting Stents Versus Durable Polymer Everolimus-Eluting Stents in Patients With ST-Segment Elevation Myocardial Infarction: 5-Year Follow-up of the BIOSTEMI Randomized Trial, Presented at: TCT 2023; October 25, 2023; San Francisco, USA. 2. Iglesias, JF. et al. Long-term outcomes with biodegradable polymer sirolimus-eluting stents versus durable polymer everolimus-eluting stents in ST-segment elevation myocardial infarction: 5-year follow-up of the BIOSTEMI randomised superiority trial, Rounded outcomes from publications, 3. Based on TLF with Orsiro DES in comparison to Xience in STEMI patient

Clinical data collected with the Orsiro DES device within the Orsiro family clinical program.

Orsiro is a trademark or registered trademark of the BIOTRONIK Group of Companies. Xience is a trademark or registered trademark of the Abbott Group of Companies.

BIOTRONIK AG  
 Ackerstrasse 6  
 8180 Bülach,  
 Switzerland

Tel +41 (0) 44 8645111  
 Fax +41 (0) 44 8645005  
 info.vi@biotronik.com  
 www.biotronik.com

© 2023 BIOTRONIK AG – All rights reserved.  
 Specifications are subject to modification, revision and improvement.

**BIOTRONIK**  
 excellence for life




# Complex Primary Percutaneous Coronary Intervention With Ultrathin Strut Biodegradable Versus Thin Strut Durable Polymer Drug-Eluting Stents In Patients With ST-segment Elevation Myocardial Infarction (STEMI): A Subgroup Analysis From The BIOSTEMI Randomized Trial<sup>1</sup>

## Conclusions

- Orsiro outperforms Xience for the treatment of non-complex primary PCI STEMI patients at 2-year follow-up (Orsiro: 4.4% ; Xience: 8.2%, p=0,008) and numerically better outcomes in complex primary PCI STEMI patients.<sup>1</sup>
- Overall, as per investigator’s interpretation, Orsiro was superior to Xience for stent-related outcomes at 2 years in primary PCI STEMI patients, irrespective of complexity.<sup>1</sup>

## Study design

Subgroup analysis of complex vs. non-complex primary PCI Acute Coronary Syndrome (ACS) patients. From the BIOSTEMI trial, an investigator-initiated, prospective, multicenter, assessor-blinded, randomized (1:1), controlled, superiority trial comparing Orsiro and Xience in STEMI patients undergoing primary PCI.




## Endpoints

### Primary Endpoint

Target Lesion Failure (TLF) at 24-month follow-up defined as the composite of:

- Cardiac Death
- Target Vessel-Myocardial Re-Infarction (TV-reMI)
- Clinically Indicated-Target Lesion Revascularization (CI-TLR)

## Characteristics of complex primary PCI patients




## Baseline clinical characteristics

Baseline clinical characteristics	Complex primary PCI		Non-complex primary PCI	
	Orsiro n = 200	Xience n = 221	Orsiro n = 653	Xience n = 626
Patients				
Age, years*	63.4 ± 11.7	64.4 ± 12.0	61.5 ± 12.0	62.3 ± 12.0
Male gender	80.5%	74.2%	78.9%	74.1%
Body mass index, kg/m <sup>2</sup> *	27.1 ± 4.7	27.3 ± 4.1	26.9 ± 4.2	26.7 ± 4.3
Diabetes mellitus	15.6%	11.8%	10.9%	13.3%
Hypertension	49.0%	51.8%	43.2%	45.0%
Hypercholesterolemia	55.3%	50.5%	46.8%	47.1%
Current smoker	44.1%	37.9%	46.6%	39.9%
Previous MI	2.5%	3.6%	4.6%	4.0%
Previous PCI	3.5%	4.1%	4.9%	5.3%
Previous CABG	1.5%	0.5%	0.6%	1.3%
Atrial fibrillation	4.0%	2.7%	1.5%	2.7%
Previous stroke or TIA	2.0%	1.8%	2.1%	2.9%
Left ventricular ejection fraction*	47.6 ± 11.3 <sup>a</sup>	46.9 ± 11.0 <sup>b</sup>	49.6 ± 10.8 <sup>c</sup>	49.0 ± 11.2 <sup>d</sup>
Multivessel disease	63.3% <sup>e</sup>	63.1% <sup>f</sup>	45.4% <sup>g</sup>	43.4% <sup>h</sup>


\* Data shown as mean ± SD; <sup>a</sup> n = 155; <sup>b</sup> n = 180; <sup>c</sup> n = 403; <sup>d</sup> n = 381; <sup>e</sup> n = 147; <sup>f</sup> n = 160; <sup>g</sup> n = 496; <sup>h</sup> n = 451

# Clinical outcomes at 2 years<sup>1</sup>


## TLF




## Cardiac Death



## TV-MI



## CI-TLR



Orsiro outperforms Xience for the treatment of non-complex, and shows numerically better outcomes in complex primary PCI STEMI patients at 2 years.

	Complex primary PCI				Non-complex primary PCI				
	Orsiro n = 200	Xience n = 221	HR [95% CI]	p-value	Orsiro n = 653	Xience n = 626	HR [95% CI]	p-value	p-value for interaction
<b>Target Lesion Failure<sup>a</sup></b>	14 (7.1%)	25 (11.6%)	0.62 (0.32-1.19)	0.15	28 (4.4%)	49 (8.2%)	0.54 (0.34-0.86)	0.008	0.74
Cardiac Death	8 (4.0%)	8 (3.6%)	1.11 (0.41-2.95)	0.84	14 (2.2%)	22 (3.6%)	0.61 (0.31-1.19)	0.14	0.32
MI	5 (2.7%)	13 (6.2%)	0.42 (0.15-1.19)	0.09	23 (2.7%)	17 (2.9%)	1.29 (0.69-2.42)	0.43	0.06
TV-reMI	2 (1.0%)	9 (4.3%)	0.24 (0.05-1.13)	0.051	10 (1.6%)	10 (1.7%)	0.95 (0.39-2.28)	0.90	0.12
Cardiac death or any MI	12 (6.2%)	21 (9.6%)	0.63 (0.31-1.28)	0.20	36 (5.7%)	38 (6.3%)	0.90 (0.57-1.43)	0.66	0.40
<b>Any Revascularization</b>	11 (5.9%)	30 (14.3%)	0.39 (0.20-0.79)	0.006	40 (6.5%)	44 (7.5%)	0.86 (0.56-1.32)	0.50	0.056
Any TLR	7 (3.7%)	17 (8.2%)	0.45 (0.19-1.10)	0.07	17 (2.8%)	26 (4.5%)	0.62 (0.33-1.14)	0.12	0.58
CI-TLR	7 (3.7%)	16 (7.7%)	0.48 (0.20-1.18)	0.10	15 (2.4%)	25 (4.3%)	0.57 (0.30-1.07)	0.08	0.78
Any TVR	9 (4.8%)	24 (11.5%)	0.41 (0.19-0.88)	0.02	23 (3.7%)	30 (5.2%)	0.72 (0.42-1.25)	0.24	0.23
CI-TVTR	9 (4.8%)	23 (11.0%)	0.43 (0.20-0.92)	0.03	21 (3.4%)	29 (5.0%)	0.68 (0.39-1.20)	0.18	0.33
<b>Target Vessel Failure<sup>b</sup></b>	16 (8.2%)	32 (14.8%)	0.54 (0.30-0.99)	0.043	36 (5.7%)	54 (9.0%)	0.63 (0.41-0.96)	0.03	0.70
POCE <sup>c</sup>	25 (12.7%)	40 (18.3%)	0.67 (0.41-1.11)	0.12	64 (10.1%)	72 (11.8%)	0.86 (0.61-1.20)	0.36	0.44
Cerebrovascular event (any)	2 (1.1%)	5 (2.3%)	0.44 (0.09-2.28)	0.32	8 (1.3%)	10 (1.7%)	0.76 (0.30-1.93)	0.56	0.57
Def. Stent Thrombosis	2 (1.0%)	6 (2.8%)	0.37 (0.07-1.84)	0.21	8 (1.3%)	8 (1.3%)	0.95 (0.36-2.53)	0.92	0.32
Def./ Prob. Stent Thrombosis	4 (2.0%)	12 (5.5%)	0.37 (0.12-1.14)	0.07	12 (1.9%)	13 (2.2%)	0.87 (0.40-1.92)	0.74	0.21

a. Composite of cardiac death, target vessel myocardial reinfarction (Q-wave and non-Q-wave), and clinically indicated target lesion revascularization (primary endpoint); b. Composite of cardiac death, any myocardial reinfarction, or any target vessel revascularization; c. Composite of all cause death, any myocardial reinfarction, or any revascularization.

## Principal investigators

Dr. J.F. Iglesias, Geneva University Hospital, Geneva, Switzerland

Dr. O. Muller, Lausanne University Hospital, Lausanne, Switzerland

1. Iglesias et al. Complex primary percutaneous coronary intervention with ultrathin-strut biodegradable versus thin-strut durable polymer drug-eluting stents in patients with ST-segment elevation myocardial infarction: A subgroup analysis from the BIOSTEMI randomized trial, Catheter Cardiovasc. Interv., 2023.

Clinical data collected with the Orsiro DES device within the Orsiro family clinical program.

Orsiro is a trademark or registered trademark of the BIOTRONIK Group of Companies. Xience is a trademark or registered trademark of the Abbott Group of Companies.

BIOTRONIK AG  
Ackerstrasse 6  
8180 Bülach,  
Switzerland

Tel +41 (0) 44 8645111  
Fax +41 (0) 44 8645005  
info.vi@biotronic.com  
www.biotronic.com

© 2023 BIOTRONIK AG – All rights reserved.  
Specifications are subject to modification, revision and improvement.

 **BIOTRONIK**  
excellence for life